

TOPAZ[®] PHD

Simple and portable video magnifier,
brought to you by

Optelec



**Designed for everyday living:
reading, writing and self-viewing**

TOPAZ® PHD

The TOPAZ PHD portable video magnifier is designed for your lifestyle. Easily carry it around the home, and between locations. Simply open, unfold the screen, switch-on and you are ready to use within seconds.

Enjoy the freedom to **READ** and **WRITE**

Once opened, the TOPAZ PHD provides ample space beneath the camera for writing and practical activities. The 15-inch screen displays more words per line and magnifies materials up to 31 times.

Check your appearance at the press of a button

Simply switch to the self-viewing camera and view yourself on-screen. Adjust the zoom to see yourself more clearly when shaving or applying make-up.

Key benefits

- Continuous zoom (2.3x – 31x)
- Vibrant 15-inch widescreen
- Self-viewing camera (up to 2.5x)
- Lightweight, foldable and portable
- High contrast semi-colours for more fluent reading (33 colour options)
- Easy to see, tactile controls
- Open, unfold the screen and switch on
- Position locator, reading lines/masks
- 4-hours continuous use with the rechargeable battery
- Carrycase included



View yourself on-screen



Plenty of space beneath the camera for reading and writing

To request a home demonstration, call

Freephone

080 8090 8090

For more information and a
You **Tube** **presentation**

Visit our website

www.optelec.co.uk

Optelec

— Part of —
VFO