


NEW!



ONYX[®] OCR




Transportable 4-in-1 Magnifier with Speech



**Simple and effective solutions for the
blind and visually impaired**

Touchscreen Point & Read



*Tactile controls are
easy to see*



To request a home demonstration of the O

Do you find reading tiresome? Would you like to combine magnification and speech in a single solution?

Now you can with the ONYX® OCR – the Transportable Desktop Magnifier with Audio Playback. With this 4-in-1 solution, you can magnify and read from newspapers and magazines on your desk, view images at a distance, view yourself, and even listen to documents as they are read aloud.

Imagine focusing the camera towards your desktop and reading from a daily newspaper, completing a crossword or enjoying craftwork, and all in High Definition. Use the tilt-and-swivel flexible camera to focus at a distance or point towards yourself to check your appearance.

With powerful magnification at your fingertips, and a wide choice of high contrast semi-colour combinations, including Yellow-text-on-black, the ONYX® OCR is your all-in-one magnification and reading solution.

Enjoy listening to long documents with Text-to-Speech

When reading (using magnification) becomes too tiring, simply sit back, relax and listen as everyday printed items are scanned and read aloud, in an instant.

Simply place a magazine or letter beneath the OCR scanning camera on the left, swipe the

screen to activate the speech features, and enjoy listening in a natural sounding voice, without straining your eyes.

With **Point & Read**, you only have to swipe your finger on the touchscreen to instantly listen to printed letters, books or magazine articles - by paragraph, column or word.

The convenient carry handle provides the possibility to move from room-to-room, around the home, or 'hotdesk' around the workplace or whilst studying.



Check your appearance or view at a distance with the Flexible camera

Key benefits

- Transportable design with carry handle
- 24-inch full colour TFT widescreen
- Working height beneath camera: 41.6 cm / 16.4 inches
- 11.7 kg / 25.8 lbs



Live Magnification

- 4-in-1 autofocus High Definition camera for desktop, distance and self-viewing
- Powerful magnification: 1.7x to 131x
- High contrast colours for fluent reading
- Full colour display for viewing photographs and magazines
- Tactile, easy to see controls



Audio, Text-to-Speech

- Powerful OCR camera (13 megapixels) which automatically recognises columns
- Scans documents up to A4 size – simply swipe the 24-inch touchscreen to begin
- Starts reading aloud within seconds
- Easy touchscreen controls for play, pause, reading speed, volume and zoom
- High contrast viewing: select preferred text and background

Optelec Limited

Unit H,
Penfold Trading Estate,
Imperial Way,
Watford WD24 4YY.
Freephone: 080 8090 8090
info@optelec.co.uk
www.optelec.co.uk



Follow Optelec on...



Request a Home Demo

call Freephone 080 8090 8090