



magnification *Optelec Traveller HD*



Key benefits

- Read at a **table**, in an **armchair** or in **bed**
- Easy to **take with you** and **store away**
- **Slide & Read:**
 - **Natural** reading experience
 - **Line marker** for easier reading
- **More words on-screen** for easier recognition
- **Continuous zoom:** 3 - 30X
- **13.3"** widescreen
- **High Definition camera** for superior images
- **Snapshot** and save up to 1,000 images
- **Rechargeable** battery offering 3 hours continuous use
- Weight: **less than 2 kg**

Optelec Limited

2 Millfield House, Woodshots Meadow,
Croxley Green Business Park,
Watford WD18 8YX.

Phone: 01923 23 13 13

info@optelec.co.uk

Buy online: www.optelec.com

Optelec



magnification




speech



braille

Traveller HD



magnification

Just read

Optelec

Simple and effective solutions for the
blind and visually impaired



Optelec Traveller HD

With the Traveller HD, you can read anywhere: at a table, on the couch, in an armchair, and even in bed! It is small and portable, so you can easily take it with you to use at home, in the office, at school, on holiday, or any other place.

Take it anywhere

Unlike desktop magnifiers, the new Traveller HD allows you to comfortably read letters, newspapers and other materials at the dining table, in bed, or in your favourite chair. Wherever you are, at home, in the office, at school or on holiday, the Traveller HD is your ideal light weight reading companion.

Superior image quality

In high definition, enjoy stepless and variable zoom, with adjustable magnification starting at 3 times, the Traveller HD allows you to view text and smaller details of photos on the high brightness 13.3-inch widescreen.

Slide & Read

Improve your reading fluency with the Traveller HD. Its unique **Slide & Read** mechanism provides you with a more natural reading experience, especially for those people that find it difficult to keep their hands steady. Simply slide the screen across a line of text and view the content on-screen. It has never been easier to follow the line or sentence. Sliding back and forth enables you to recognise the shapes of individual words. You can see a more complete line of text because of the large screen, so you can read more fluently.

