

# ClearView C Speech

X-Y reading table for fluent reading

Height-adjustable tilting screen



*Sit back and relax as printed documents are magnified and read aloud.*



**Simple and effective solutions for the  
blind and visually impaired**



## Home demo? Call 080 8090 8090

The **ClearView C Speech** makes reading easier. It combines a clear voice with magnification and high contrast. When reading becomes too tiring, simply activate the speech feature and enjoy listening to magazine articles, correspondence and books.

### Key benefits

- Magnify letters and magazines from 1.5x - 75x or activate the Speech Features (listed below)
- Scans documents up to A4 size – simply swipe the 24-inch touchscreen to begin
- Starts reading aloud within seconds with **Full Page Overview**
- **Point & Read**: touch a column, paragraph or word on-screen and listen instantly
- Scroll through the content with your finger, even when magnified (from 1 to 75x)
- Choose from 74 voices/accents in 36 reading languages
- **Archive documents/photos**: view and listen at your convenience
- **Picture Viewer**: zoom and slideshow options to view photos
- Built-in stereo speakers and headphone connection

*Optelec Limited*  
Unit H,  
Penfold Trading Estate,  
Imperial Way,  
Watford WD24 4YY.  
Freephone: 080 8090 8090  
info@optelec.co.uk  
www.optelec.co.uk



Follow Optelec on...

